



Roscoe Economy Knee Scooter

The **Roscoe Economy Knee Scooter** offers an affordable choice in mobility for those suffering from an ankle or foot injury. The Roscoe Economy Knee Scooter can easily be folded down for transportation and storage. The handle bars and knee platform are height adjustable to accommodate a wide range of users. The frame and knee pad are manufactured for worry free use.



Height Adjustable Steering Wheel



Adjustable Hand Brake

FEATURES

- Scooter can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Quick Release Steering Assembly folds for easy storage
- Quick Release Seat includes height adjustments to accommodate different patient heights
- Padded Knee Platform with easy adjustment mechanism
- Oversize wheels are ideal for indoor and outdoor use
- Adjustable braking system
- Limited lifetime warranty on frame, 90 days on Non-Durable components

COLORS

Sterling Grey



SPECIFICATIONS

Patient Height Range:	4'11" – 6'6"
Adjustable Handle Height:	31" – 41"
Product Weight:	19.5 lbs.
Product Dimensions:	31.75" x 16.125" x 30"
Weight Capacity:	250 lbs.

Item #	Description	UPC	UOM
ROS-KS2	Roscoe Economy Knee Scooter	092237614624	each