

Electrotherapy InTENSity™ Select Combo

The InTENSity™ Select Combo is a powerful, easy-to-use portable device. Complete with four waveforms, the unit gives you and your patients the flexibility to treat a variety of musculoskeletal issues. With a strong current of 105 mA, this innovative device combines TENS, EMS, Interferential and Microcurrent and comes complete with an AC adaptor ensuring you always have power to the device when traveling.



Product Features:

- Simple to operate, portable device for pain management
- Four waveforms, TENS, EMS, Interferential and Microcurrent
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, AC adaptor, manual and pack of four self-adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adaptor
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb. with battery
- One year warranty
- HCPCS Code: E0730* for TENS Unit
- HCPCS Code: A4595* for electrotherapy supplies



Product Specifications:

	TENS	EMS	IF	MICRO
Wave Form	Monophasic Square	Biphasic	Biphasic Square	Monophasic Square
Pulse Rate	1-150 Hz	1-150 Hz	4kHz	1-150 Hz
Pulse Width	50-300 μs	50-300 μs	125 μs, phase duration 250 μs	2-200 μs
Output Voltage	Constant Current	Constant Current	Constant Current	Constant Current
Current	0-105 mA	0-105 mA	0-105 mA	0-700 μA
Modes	4	3	4	3
Timer	1-60 Constant	1-60 Constant	1-60 Constant Interf. Frequency 1-150 beats	1-60 Constant

Item #	Description	Qty
D18195	InTENSity™ Select Combo	each

* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com